

Specialized Injector for Body, Scalp and Face

Fine, Clear Injection amount controling Controling Range: 0.01cc ~ 0.2cc / Continue 1 Setting, Same amount injection



Features

- · Specialized for Body, Face and Scalp Solution Injection
- · Drug Loss Control (Patent)
- · Air Compressor type : Powerful Injection, Minimize the pain
- · Available for Body, Face and Scalp Meso theraphy

Lipolysis program



step 1

RelaxationRF system

step 2

Lipolysis mesotherapy

Airmeso

step 3

Blood circulation

LED therapy

step 4

Oxidation

Combination physical therapy

Scalp therapy program



step 1

Blood circulation

Hair steamer

step 2

Peeling & Deep cleansing

Peeling machine

step 3

Scalp nutritious mesotherapy

Airmeso

step 4

Soothing, Blood circulation


LED therapy



BEST DRUG INJECTION EQUIPMENT!

Popularization of drug injector faithful to performance and reasonable price.









Mode & Display

0.01cc [0]

Dose mode

Inject set amount solution by shot 0.01cc ~ 0.2cc



Nappage mode

Automatically inject set amount solution quickly 0.01cc ~ 0.1cc (10 Level)



Auto mode

Fast needling for stimulating 100 ~ 500 shot/min (5 Level)



Continue mode

Inject solution continuously

Specification

- Rated Voltage: AC 220V/50~60Hz/90W
- Gun size: 264mm X 282mm X 243mm
- · Gun weight: 600g
- Body size: 300mm X 350mm X 250mm
- · Body weight (including Gun): 9.0kg
- Mode: Dose mode (0.01~0.2cc)
 Nappage mode (0.01~0.1cc)
 Auto mode (100~500 shot/min)
 Continue mode

